

Overview of CoolBox XT Workstation

CoolBox XT is a portable workstation designed for maintaining sample temperatures below 4°C on the bench top without the use of wet ice or electricity. Patent-pending dual-phase conductive XT Cooling Core or XT Freezing Core provides the cooling source when a thermo-conductive CoolRack or CoolSink tube or plate module is placed on top. The core and sample module in combination ensure uniform well-to-well temperature throughout the cooling period regardless of sample position. XT Cooling Core can maintain sample temperature from 0.5 to 4°C for over 16 hours; XT Freezing Core can maintain frozen samples for up to 10 hours; using dry ice in the base provides an ultra-cold temperature (-78°C) for snap-freezing samples in tubes or plates. For a list of CoolRack and CoolSink modules that are compatible with the CoolBox XT, see below.

Working Temperature Range	CoolBox XT Cooling Source	Temperature Duration* Open Lid	Temperature Duration* Closed Lid
0.5 to 4°C	XT Cooling Core	Over 10 hours	Over 16 hours
-20 to 0°C	XT Freezing Core	Over 5 hours	Over 8 hours
-78°C	200 ml Dry Ice	Over 4 hours	Over 5 hours

*All tests were performed using a CoolRack XT M24 loaded with 24 TruCool® 2.0ml microcentrifuge tubes filled with 1.5 mL water. Actual performance may vary depending upon CoolRack module employed, sample load, initial sample temperature, ambient temperature, air currents, and other conditions.

Quick Start

- Remove XT Cooling Core from -20°C freezer and place on bench
- When temperature strip on XT Cooling Core registers 1°C, place CoolRack® or CoolSink® sample module in the CoolBox XT base and place collar on top
- Load samples in CoolRack or CoolSink modules

Caution

- Lift CoolBox XT from the bottom using the handholds on the sides.
- Always use two hands when carrying or lifting CoolBox XT.
- Avoid touching the top metal surface of the XT Cooling Core and XT Freezing Core when removing from freezer.


IMPORTANT: To ensure your samples do not freeze and to get maximum cooling duration, please refer to the detailed instructions that follow.

Assembly of CoolBox XT

The CoolBox XT System consists of a high-density polyethylene foam base, collar and lid, and a reusable XT Cooling Core. A CoolBox XT Workstation also includes a CoolRack or CoolSink tube or plate module. When placed on top of the XT Cooling or Freezing Core, the modules and samples rapidly equilibrate to the temperature of XT Core. Samples can be placed in the CoolRack or CoolSink module either before or after the module is placed in CoolBox XT. When used with CoolSink plate modules the magnetized collar is not required and can be removed for easy access when working; however, to obtain the maximum cooling duration, we recommend using the collar whenever possible.

Assembly of CoolBox XT cont.

1. Remove XT Cooling Core from the freezer and place on benchtop for approximately 10-15 minutes. When the temperature indicator displays 1°C, XT Cooling Core is ready to use. *Note: Frost will form on the core exterior upon removal from the freezer; when the frost liquifies, the core is at proper temperature*
2. Place XT Cooling Core into CoolBox XT base.
3. Fit magnetized collar onto CoolBox XT base, seating it securely.
4. Place the CoolRack or CoolSink module of choice onto XT Cooling Core.
5. Load samples.
6. Place lid on CoolBox XT when not processing samples to maximize cooling duration.



Using XT Cooling Core for Maintaining Samples at 0.5 to 4°C

Working Temperature Range	CoolBox XT Cooling Source	Temperature Duration Open Lid	Temperature Duration Closed Lid
0.5 to 4.0°C	XT Cooling Core	Over 10 hours	Over 16 hours

Freeze XT Cooling Core in a -20°C freezer for at least 12 hours. XT Cooling Core should be stored in -20°C freezer when not in use so it is ready when needed. *Note: Freezing XT Cooling Core for less than the specified time will result in decreased cooling duration.*

When using a CoolBox XT to cool samples (1 mL or greater) in a CoolRack using a charged XT Cooling Core

- Remove XT Cooling Core from the freezer and place into base.
- Fit magnetized collar onto base, seating it securely.
- Place CoolRack or CoolSink module directly onto core and allow to equilibrate to 4°C (approximately 10 -15 minutes).
- Load samples
- Place lid on CoolBox XT when not processing samples to maximize cooling duration.

When using a pre-chilled (4°C) CoolRack or CoolSink module

- Remove XT Cooling Core from the freezer and place on benchtop for approximately 10 minutes. When the temperature indicator displays 1°C, XT Cooling Core is ready to use.

IMPORTANT: Failure to allow XT Cooling Core to reach 1°C may result in sample freezing.

- Place XT Cooling Core into base.
- Fit magnetized collar onto base, seating it securely.
- Place the CoolRack or CoolSink module of choice onto the core.
- Load samples.
- Place lid on CoolBox XT when not processing samples to maximize cooling duration.

Using Optional XT Freezing Core for Maintaining Samples at -20 to 0°C

Working Temperature Range	CoolBox XT Cooling Source	Freezer Temperature	Temperature Duration* Open Lid	Temperature Duration* Closed Lid
-20 to 0°C	XT Freezing Core	-20°C	Over 5 hours	Over 8 hours
-20 to 0°C	XT Freezing Core	-80°C	Over 8 hours	Over 12 hours

Freeze XT Freezing Core in a -20°C freezer for at least 12 hours, or for a faster start, freeze in a -80°C freezer for at least 6 hours. Freezing in a -80°C freezer also prolongs cooling duration. XT Freezing Core should be stored in a -20°C or -80°C freezer when not in use so it is ready when needed.

Note: Freezing XT Freezing Core for less than the specified time will result in decreased cooling duration.

When using a room temperature CoolRack or CoolSink module

- Remove XT Freezing Core from the freezer and place into base.
- Fit magnetized collar onto the base, seating it securely.
- Place the room temperature CoolRack or CoolSink module of choice onto core and allow to equilibrate to 0°C (approximately 10 minutes).
- Load samples.
- Place lid on CoolBox XT when not processing samples to maximize cooling duration.

When using a pre-chilled (10°C or less) CoolRack or CoolSink module

- Remove XT Freezing Core from the freezer and place into base.
- Fit magnetized collar onto the base, seating it securely.
- Place the CoolRack or CoolSink module of choice onto core.
- Load samples.
- Place lid on CoolBox XT when not processing samples to maximize cooling duration.

Using Dry Ice for Maintaining or Snap-Freezing Samples at -78°C

Working Temperature Range	CoolBox XT Cooling Source	Temperature Duration* Open Lid	Temperature Duration* Closed Lid
-78°C	200ml dry ice	Over 4 hours	Over 5 hours

- Remove XT Cooling or Freezing Core from the CoolBox XT base.
- Fill the base, with approximately 200ml of pulverized dry ice.
- Fit magnetized collar onto the base, seating it securely.
- Place the CoolRack or CoolSink module directly onto dry ice and allow module to equilibrate to -78°C (approximately 7-8 minutes).
- Load samples.
- If snap-freezing, freezing will occur in 7-10 minutes depending upon sample volume and type.

Note: the thermo-conductive design of the CoolRack and CoolSink modules ensures uniform well-to-well temperature regardless of the consistency of the dry ice.

CoolBox XT Care and Cleaning

The CoolBox XT housing is constructed from a cross-linked closed-cell dense polyethylene foam. The material has excellent resistance to fluid absorption and abrasion. Do not use the CoolBox XT base for pulverizing dry ice. Maximum temperature exposure: 60°C. Avoid prolonged exposure to UV light sources.

All components including housing, XT Cooling Core and optional XT Freezing Core are compatible with repeated and prolonged cryogenic temperature exposure. All components can be cleaned with aqueous detergents, alcohol, 10% bleach, and acid/base viricide (such as Virkon S) solutions. Rinse with clear water after using cleaning solutions. Do not autoclave.

CoolRack and CoolSink sample modules may be autoclaved, or cleaned with alcohol or 10% bleach.

CoolBox XT Dimensions

Interior (with XT Cooling Core in base): 14 x 9.4 x 9.9 cm (L x W x H) / 5.5 x 3.7 x 3.9 in (L x W x H)

Exterior: 20 x 16 x 15 cm (L x W x H) / 7.8 x 6.3 x 5.9 in (L x W x H)

CoolRack and CoolSink Thermo-Conductive Sample Modules Compatible with CoolBox XT

Item	Description	Accommodates
BCS-163	CoolRack M6	6 x 1.5ml or 2ml microfuge tubes
BCS-125	CoolRack M15	15 x 1.5ml or 2ml microfuge tubes
BCS-127	CoolRack M15-PF	15 x 1.5ml conical microfuge tubes
BCS-535	CoolRack XT M24	24 x 1.5 or 2.0 mL microcentrifuge tubes, SBS-compliant
BCS-126	CoolRack CF15	15 x 1ml or 2ml cryovials
BCS-534	CoolRack XT CFT24	24 cryogenic vials, with "gripping" wells for one-hand vial open/close, SBS-compliant
BCS-529	CoolRack XT PCR96	PCR tubes, strips or 96-well plate
BCS-523	CoolRack XT M-PCR	12 x 1.5 or 2.0 mL microcentrifuge tubes and 6 PCR strip wells
BCS-538	CoolRack XT PCR384	one 384-well PCR plate, SBS-compliant
BCS-231	CoolRack 96x0.5ml	96 x 0.5ml 2D barcode tubes
BCS-149	CoolRack 96x1ml	96 x 1ml 2D barcode tubes
BCS-536	CoolSink XT 96F	one flat-bottom 96- or 384-well assay plate, SBS-compliant
BCS-537	CoolSink XT 96U	one u-bottom 96-well assay plate, SBS-compliant
Tall Tube Modules - requires Extension Collar, if lid closing is necessary		
BCS-232	CoolRack L	12 x 15 mL centrifuge tubes
BCS-153	CoolRack 15ml	9 x 15 mL centrifuge tubes
BCS-154	CoolRack 50ml	4 x 50 mL centrifuge tubes
BCS-235	CoolRack LV	12 x 13mm or 16mm diameter tubes
BCS-157	CoolRack VS13	9 x 13x75mm blood tubes
BCS-155	CoolRack V13	9 x 13x100mm blood tubes or 5ml cryogenic vials
BCS-156	CoolRack VS16	9 x thermo-conductive module for 16x100mm blood tubes

⚠ CAUTION: *The products described here are intended for the exclusive use by trained and experienced laboratory and medical personnel. Use of dry ice can be dangerous. Direct skin contact with dry ice or metal components that have been in contact with dry ice can cause freezing injury. Always use appropriate protective equipment for eyes and skin when handling dry ice and cold metal components.*

Made in USA.

biocision®

12 E. Sir Francis Drake Blvd, Suite B, Larkspur CA, 94939 USA info@biocision.com www.biocision.com

CoolBox XT Instructions, January 2014
CoolBox XT ©2014. BioCision, LLC. All rights reserved. Patents pending. BioCision, Standardizing Samples CoolCell, CoolBox, CoolRack, CoolSink designations are trademarks owned by BioCision, LLC.