

Ag+goods!

Silver Ion Water Preservation Cell

Ag+goods! is a silver anti-microbial agent. Just put it into the water, Ag+goods! keeps water clean and stops growth of microorganisms.

Ag+goods! is suitable for CO2 incubators, Water baths and Water chillers.

Details:

■ Easy to Use

Keeps water clean, just put Ag+goods! into water. It remains effective for up to six months.

■ Power of Silver Ion(Ag+)

Ag+goods! reduces the risks of contamination by bacteria, viruses and other microorganisms.

■ Special Structure

Even if water level drops, Ag+goods! remains effect because of its thin designing.

■ Harmless for Cells


The silver ions do not diffuse into the air. Ag+goods! keeps water clean and has no effect on cultured cells.

■ High Antimicrobial Effect

Performance tests carried out with *Staphylococcus aureus*, *E coli* and viruses showed that Ag+goods! kills them within 8 hours of its being placed in contaminated water.



Example use in a CO₂ incubator



Performance test with *Staphylococcus aureus*

WAKEN  TECH CO., LTD.

WAKENBTECH CO.,LTD.

1-12-8, Senbahigashi, Minoh-city, Osaka 562-0035, JAPAN

T: +81-72-749-5300 F: +81-72-749-5600

E: info@wakenbtech.co.jp W: www.wakenbtech.co.jp