

Edgeless Rack CRA-5006

Tube rack for 6x50mL centrifuge tubes

This edge-less rack with a unique design keeps tubes safe and secure when they are being handled in the laboratory.

The rack holds six 50mL tubes and the absence of sharp edges prevents accidents and makes the rack very easy to clean.

Details:

■ Unique Design

This product is made of uncoated stainless steel wire. It has no sharp edges or corners, is easy to load and safe to use.

■ No scratching

No sharp edges to scratch the hands, or to tear or puncture laboratory gloves.

■ Clean Design

Proper intervals between tubes prevent cross-contamination.

■ High Usability

Tubes are easy to place and remove.


■ Design does not allow fibers to accumulate

This rack is easy to clean and will not accumulate fibers from laboratory wipes.

■ Suitable for Clean rooms and Isolators

The clean and safe design is suitable for use in an aseptic environment.

Problems of conventional tube racks



Burrs



Fibers Accumulation



Tube caps conflict next ones


Advantage of Edgeless Rack



Design without burr.



Proper intervals between tubes.



Easy to get out / in tubes.

WAKEN  TECH CO., LTD.

WAKENBTECH CO.,LTD.

1-12-8, Senbahigashi, Minoh-city, Osaka 562-0035, JAPAN

T: +81-72-749-5300 F: +81-72-749-5600

E: info@wakenbtech.co.jp W: www.wakenbtech.co.jp