


ELISA / Western blotting blocking agent

Improve your results by eliminating backgrounds

Block Ace

Block Ace is a new blocking agent prepared from lactprotein. Since lactprotein contains a lot of hydrophobic groups , it shows a high affinity to plastic , membrane filter , etc . Block Ace exhibits an more excellent blocking effect on nonspecific adsorption of antigens or antibodies in the ELISA , RIA and immunoblocking methods compared with the conventional blocking agents such as BSA(Bovine Serum Albumin)

- **Easy to use**
- **Reduce non-specific binding**
- **Reliable consistent formulation**



Superior to 1% BSA , BLOCK ACE is a high-performance blocking reagent for ELISA assays and Western blotting . In comparison to BSA , BLOCK ACE provides reduced backgrounds and sharper standard curves.It can also be used as a wash buffer , and for diluting antibodies .

Description	Quantity	Product Code
BLOCK ACE	20×4g	UK-B80
	1 ×40g	UK-B40
	1 ×500g	UK-B500

【Manufacturer】



KAC Co.,LTD

【Distributed by】